

Saving your life is a way of their lives

The Saviour Medical Rescue Stretcher is the latest versatile rescue stretcher to be available to the emergency, resilience and rescue services as a single product solution. Ideally suited for work in confined spaces, rescue, hazardous and water based environments, the stretcher is an innovative design likened to being the “Swiss Army Knife” of rescue stretchers.

This lightweight device (5kg) with its polyethylene skin and removable straps is infection control compliant. The stretcher is neutrally buoyant making water based rescue possible and the sand and waterproof design of the fixings further support such use. The stretcher, fixings and strap system are load rated to 200kg (31 stone) making it suitable for most patients. The integral headblock and double skin in the spinal area assists in the management of patients with potential spinal injuries.

The stretcher is supplied as a rolled device in a robust carry bag, making it compact and aiding access to confined spaces as well as requiring little storage space. The device is also X-Ray translucent and CT friendly. The simplistic, colour coordinated strap system ensure there is little additional training requirement. The stretcher is CE marked and manufactured in the UK.

The Saviour Medical Rescue Stretcher has already been rigorously tested across the UK by a renowned water based rescue organisation. The tactical variant of the stretcher is also in use by the military.

The UK emergency services, local Government and national rescue organisations have experienced increased demand for suitable rescue equipment to ensure emergency preparedness and resilience.

The Saviour Medical Rescue Stretcher has been designed and developed by a group of highly skilled and experienced rescuers. The stretcher is manufactured in the UK and distributed through DS Medical, a Hampshire based medical supplier.

Contact

Mr Richard Bourke:
07960 866335

Mr Jason Warrener:
07872 039998

Data Southern (Enterprise) Ltd
(T/A DS Medical)
Unit 6 Fareham Enterprise Centre
Hackett Way
Fareham
PO14 1TH
Tel.: **01329 311451**
Fax: **01329 829276**
Email: info@dsmedical.co.uk
Web: www.dsmedical.co.uk

